

## feline IFN-gamma/IFNG protein (Recombinant) (His-Tag) (STJP019304)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	19.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Immunogen Region	24-167aa
Immunogen Sequence	MGSSHHHHHH SSSLVPRGSH MGSQAMFFKE IEELKGYFNA SNPDVADGGS LFDILKNWK EESDKTIQS QIVSFYKMF ENLKDDDDQRI QRSMDTIKED MLDKLLNTSS SKRDDLKLI QIPVNDLQVQ RKAINELFKV MNDLSPRSNL RKRKRSQNLF RGRNRASK

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant feline IFN-gamma IFNG protein under reducing conditions stained with Coomassie blue