

Human DUSP18 protein (Recombinant) (His-Tag) (STJP019320)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 40% Glycerol, 0.1mM PMSF, 1mM EDTA
Molecular Weight	23.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	150290
Gene Symbol	DUSP18
UniProt ID	DUS18_HUMAN
Immunogen Region	1-188aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSHTAPSC AFPVQFRQPS VSGLSQITKS LYISNGVAAN NKMLLSSNQI TMVINVSVEV VNTLYEDIQY MQVPVADSPN SRLCDFFDPI ADHIHSVEMK QGRDLLHCAA GVSRSALCL AYLMKYHAMS LLDHTWTKS CRPIIRPNSG FWEQLIHYEF QLFGKNTVHM VSSPVGMIPD IYEKEVRLMI PL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



15% SDS-PAGE (3µg)

SDS-PAGE gel of recombinant human DUSP18 protein 3µg under reducing conditions stained with Coomassie blue