

## Human VAMP-2 protein (Recombinant) (His-Tag) (STJP019340)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol, 0.1mM PMSF
Molecular Weight	12.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">6844</a>
Gene Symbol	<a href="#">VAMP2</a>
UniProt ID	<a href="#">VAMP2_HUMAN</a>
Immunogen Region	1-94aa
Immunogen Sequence	MSATAA TAPPAAPAGE GGPPAPPPNL TSNRRLQQTQ AqvDEVVDIM RVNVDKVLER DQKLSELDDR ADALQAGASQ FETSAAKLKR KYWWKNLK

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human VAMP-2 protein 3  $\mu$ g under reducing conditions stained with Coomassie blue