

Human UbchH5b/UBE2D2 protein (Recombinant) (No-Tag) (STJP019355)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM MES buffer (pH 6.0) containing 50mM NaCl, 1mM DTT
Molecular Weight	16.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7322
Gene Symbol	UBE2D2
UniProt ID	UB2D2_HUMAN
Immunogen Region	1-147aa
Immunogen Sequence	MALKRIHKEL NDLARDPPAQ CSAGPVGDDM FHWQATIMGP NDSPYQGGVF FLTIHFPTDY PFKPPKVAFT TRIYHPNINS NGSICLDILR SQWSPALTIS KVLLSICSL CDPNPDDPLV PEIARIYKTD REKYNRIARE WTQKYAM

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human UbchH5b UBE2D2 protein under reducing conditions stained with Coomassie blue