

Mouse Glutaredoxin 1/GLRX1 protein (Recombinant) (His-Tag) (STJP019358)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM MES buffer (pH 5.0) containing 20% Glycerol, 0.1M NaCl
Molecular Weight	14.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	93692
Gene Symbol	Glrx
UniProt ID	GLRX1_MOUSE
Immunogen Region	1-107aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MGSMARQEFVN CKIQSGKVVV FIKPTCPYCR KTQEILSQLP FKQGLLEFVD ITATNNTSAI QDYLLQLTGA RTVPRVFIGK DCIGGCSDLI SMQQTGELMT RLKQIGALQL

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant mouse Glutaredoxin 1 protein under reducing conditions