

Human FGF acidic/FGF-1 protein (Recombinant) (No-Tag) (STJP019385)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.0) containing 100mM NaCl
Molecular Weight	15.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2246
Gene Symbol	FGF1
UniProt ID	FGF1_HUMAN
Immunogen Region	16-155aa
Immunogen Sequence	MFNLPPGNYK KPKLLYCSNG GHFLRILPDG TVDGTDRDSD QHIQLQLSAE SVGEVYIKST ETGQYLAMDT DGLLYGSQTP NEECLFLERL EENHYNTYIS KKHAEKWV FV GLKKNKNGSCKR GPRTHYGQKA ILFLPLPVSS D

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human FGF acidic FGF-1 protein under reducing conditions stained with Coomassie blue