

## Human ICAM-5 protein (Recombinant) (hIgG-His-Tag) (STJP019439)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	HEK293 Cells

### PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	113.8kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">7087</a>
Gene Symbol	<a href="#">ICAM5</a>
UniProt ID	<a href="#">ICAM5_HUMAN</a>
Immunogen Region	23-835aa
Immunogen Sequence	LFGLSAVSQE PFWADLQPRV AFVERGGSLW LNCSTNCPRP ERGGLETSLR RNGTQRGLRW LARQLVDIRE PETQPVCFFR CARRTLQARG LIRTFQRPDR VELMPLPPWQ PVGENFTLSC RVPGAGPRAS LTLTLRGAQ ELIRRSFAGE PPRARGAVLT ATVLARREDH GANFSCRAEL DLRPHGLGLF ENSSAPREL TFSLSPDAPR LAAPRLLEVQ SERPVSCT

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human ICAM-5 protein under reducing conditions stained with Coomassie blue