

Human Resistin/RETN protein (Recombinant) (His-Tag) (STJP019450)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid In 20mM Sodium citrate (pH3.0) containing 20% Glycerol
Molecular Weight	10.3kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56729
Gene Symbol	RETN
UniProt ID	RETN_HUMAN
Immunogen Region	19-108aa
Immunogen Sequence	KTLCSMEEAI NERIQEVAGS LIFRAISSIG LECQSVTSRG DLATCPRGFA VTGCTCGSAC GSWDVR AETT CHCQCAGMDW TGARCCRVQP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human Resistin RETN protein under reducing conditions with Coomassie blue stain