

Mouse CD39L2/ENTPD6 protein (Recombinant) (His-Tag) (STJP019452)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid In Phosphate-Buffered Saline (pH 7.4) containing 20% Glycerol
Molecular Weight	47.3kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12497
Gene Symbol	Entpd6
UniProt ID	ENTP6_MOUSE
Immunogen Region	33-455aa
Immunogen Sequence	KWHRAS AAQAFFTIAG AASGARWTQQ AFSSPGSAAR GHEVFGIMF DAGSTGTRIH VFQFARPPGE TPTLTHTETFK ALKPGLSAYA DDVEKSAQGI QELLNVAKQH IPYDFWKATP LVLKATAGLR LLPGEKAQKL LQKVKEVFKA SPFLVGDDCV SIMNGTDEGV SAWITVNFLT GSLKTPGSSS VGMLDLGGGS TQITFLPRVE GTLQASPPGH LTALQM

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant mouse CD39L2/ENTPD6 protein under reducing conditions stained with Coomassie blue