

Human FOLR2 protein (Recombinant) (His-Tag) (STJP019453)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Bioactivity
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	25.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	2350
Gene Symbol	FOLR2
UniProt ID	FOLR2_HUMAN
Immunogen Region	22-230aa
Immunogen Sequence	QDRTDLLNVC MDAKHHKTKP GPEDKLHDQC SPWKKNACCT ASTSQELHKD TSRLYNFNWD HCGKMEPACK RHFIQDTCLY ECSPNLGPWI QQVNQSWRKE RFLDVPLCKE DCQRWWEDCH TSHTCKSNWH RGWDWTSGVN KCPAGALCRT FESYFPTPAA LCEGLWSHSY KVSNYSRGSG RCIQMWFD SA QGNPNEEVAR FYAAAMHVN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human FOLR2 protein SDS-PAGE gel under reducing conditions stained with Coomassie blue



FOLR2 protein binding assay with Folic Acid-BSA showing ED50 dose response curve