

Mouse TIM-3/HAVCR2 protein (Recombinant) (hIgG-His-Tag) (STJP019455)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid In Phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	46.3kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	171285
Gene Symbol	Havcr2
UniProt ID	HAVR2_MOUSE
Immunogen Region	20-193aa
Immunogen Sequence	RSLENAYVFE VGKNAYLPCS YTLSTPGALV PMCWKGKFCP WSQCTNELLR TDERNVYQK SSRYQLKGD NKGDVSLIK NVTLDDHGTY CCRIQFPGLM NDKKLELKLD IKAQVTPAQ TAHGDSTTAS PRTLTTERRG SETQTLVTLH NNGTKISTW ADEIKDSGET IRTA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant mouse TIM-3 HAVCR2 protein 3µg under reducing conditions with Coomassie blue staining