

Human C-Reactive/CRP protein (Recombinant) (His-Tag) (STJP019456)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid In Phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	24.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1401
Gene Symbol	CRP
UniProt ID	CRP_HUMAN
Immunogen Region	19-224aa
Immunogen Sequence	QTDMSRK AFVFPKESDT SYVSLKAPLT KPLKAFTVCL HFYTELSSTR GYSIFS YATK RQDNEILIFW SKDIGYSFTV GGSEILFEVP EVTVAPVHIC TSWESASGIV EFWVDGKPRV RKSLKKGYTV GAEASIILGQ EQDSFGGNFE GSQSLVGDIG NVNMWDFVLS PDEINTIYLG GPFSPNVLNW RALKYEVQGE VFTKPQLWP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human C-Reactive Protein CRP under reducing conditions
Coomassie blue stain