

Human PSGL-1/CD162 protein (Recombinant) (hlgG-His-Tag) (STJP019458)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid In Phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	53.4kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6404
Gene Symbol	SELPLG
UniProt ID	SELPL_HUMAN
Immunogen Region	42-295aa
Immunogen Sequence	QATEY EY LDYDFLPETE PPEMLRNSTD TTPLTGPGTP ESTTVEPAAR RSTGLDAGGA VTELTTELAN MGNLSTDSAA MEIQTTPAA TEAQTTPPLAA TEAQTTRLTA TEAQTTPPLAA TEAQTTPPAA TEAQTTPPTG LEAQTTPAA MEAQTTPAA MEAQTTPAA MEAQTTPAA MEAQTTPAA MEAQTTPAA TEAQTTPPTA TEAQTTPPLA MEALSTEPSA TEALSM

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human PSGL-1
CD162 protein under reducing conditions
Coomassie blue stain