

NKMAX Marker (7 bands, 7.2-70 kDa) protein (Recombinant) (STJP019468)

GENERAL INFORMATION


Product Type	Proteins
Applications	SDS-PAGE

PRODUCT PROPERTIES

Formulation	Liquid in 12mM Tris-HCl buffer (pH 6.8) 5% Glycerol, 0.4% SDS 2.88mM 2-mercaptoethanol, 0.02% bromophenol blue
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



NKMAX Protein Marker 7 bands SDS-PAGE gel under reducing conditions stained with Coomassie blue