

Human Thrombospondin 5/COMP protein (Recombinant) (His-Tag) (STJP019477)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid In Phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	82.4kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1311
Gene Symbol	COMP
UniProt ID	COMP_HUMAN
Immunogen Region	21-757aa
Immunogen Sequence	QGGQSPL GSDLGPQMLR ELQETNAALQ DVRELLRQQV REITFLKNTV MECDACGMQQ SVRTGLPSVR PLLHCAPGFC FPGVACIQTE SGARCGPCPA GFTGNGSHCT DVNECNAHPC FPRVRCINTS PGFRCEACPP GYSGPTHQGV GLAFKANKQ VCTDINECET GQHNCVPNSV CINTRGSFQC GPCQPGFVGD QASGCQRRRAQ RFCPDGSPSE CHEHAD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human Thrombospondin 5 COMP protein under reducing conditions