

Human LILRB2/CD85d/ILT4 protein (Recombinant) (His-Tag) (STJP019479)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid In Phosphate-Buffered Saline (pH 7.4) containing 20% Glycerol
Molecular Weight	48.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10288
Gene Symbol	LILRB2
UniProt ID	LIRB2_HUMAN
Immunogen Region	24-461aa
Immunogen Sequence	GTIPKPTLWA EPDSVITQGS PVTLSQCGSL EAQEYRLYRE KKSASWITRI RPELVKNGQF HIPSITWEHT GRYGCQYYSR ARWSELSDPL VLVMTGAYPK PTLAQPSPV VTSGGRVTLQ CESQVAFGGF ILCKEGEDEH PQCLNSQPHA RGSSRAIFSV GPVSPNRRWS HRCYGYDLNS PYVWSSPSDL LELLVPGVSK KPSLSVQPGP VMAPGESLTL QCVSDVGY

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant human LILRB2 CD85d ILT4 protein
 SDS-PAGE gel under reducing conditions