

E.coli rnpA protein (Recombinant) (No-Tag) (STJP019487)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid Phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	13.7kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	1-119aa
Immunogen Sequence	MVKLAFPREL RLLTPSQFTF VFQQPQRAGT PQITILGRIN SLGHPRIGLT VAKKNVRRRAH ERNRIKRLTR ESFRLRQHEL PAMDFVVAK KGVADLDNRA LSEALEKLWR RHCRLARGS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant E.coli rnpA protein 3µg under reducing conditions stained with Coomassie blue