

Human SIRP alpha/CD172a protein (Recombinant) (His-Tag) (STJP019500)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Bioactivity
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid In Phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	39kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	140885
Gene Symbol	SIRPA
UniProt ID	SHPS1_HUMAN
Immunogen Region	27-373aa
Immunogen Sequence	GVAGEEE LQVIQPKSV LVAAGETATL RCTATSLIPV GPIQWFRGAG PGRELIYNQK EGHFPRVTTV SDLTKRNNMD FSIRIGNITP ADAGTYCYVK FRKGSPDDVE FKSGAGTELS VRAPKSPV SVSGAARATPQ HTVSFTCESH GFSPRDITLK WFKNGNELSD FQTNVDPVGE SVSYSIHSTA KVLVTREDVH SQVICEVAHV TLQGDPLRGT ANLSET

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human SIRP alpha CD172a protein
 SDS-PAGE gel 3µg reducing conditions
 Coomassie stain



SIRP alpha CD172a protein binding to human
 CD47 in functional ELISA assay