

Human CENPM protein (Recombinant) (His-Tag) (STJP019502)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 30% Glycerol, 1mM DTT
Molecular Weight	22.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79019
Gene Symbol	CENPM
UniProt ID	CENPM_HUMAN
Immunogen Region	1-180aa
Immunogen Sequence	MSVLRPL DKLPGLNTAT ILLVGTEDAL LQQLADSMK EDCASELKVH LAKSLPLPSS VNRPRIDLIV FVVNLHISKYS LQNTEESLRH VDASFFLGKV CFLATGAGRE SHCSIHRHTV VKLAHTYQSP LLYCDLEVEG FRATMAQRLV RVLQICAGHV PGVSALNLLS LLRSSEGPSL EDL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human CENPM protein SDS-PAGE gel 3µg reducing conditions Coomassie stain