

Human PKC iota protein (Recombinant) (His-Tag) (STJP019507)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid In Phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	69kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5584
Gene Symbol	PRKCI
UniProt ID	KPCI_HUMAN
Immunogen Region	1-596aa
Immunogen Sequence	MPTQRDSSTM SHTVAGGGSG DSHQVRVKA YYRGDIMITH FEPSISFEGL CNEVRDMCSF DNEQLFTMKW IDEEGDPCTV SSQLELEEF RLYELNKDSE LLIHVFPVCP ERPGMPCPGE DKSIYRRGAR RWRKLYCANG HTFQAKRFNR RAHCAICTDR IWGLGRQGYK CINCKLLVHK KCHKLVTEIC GRHSLPQEPV MPMDQSSMHS DHAQTVIPYN PSSHESLD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human PKC iota protein under reducing conditions with Coomassie stain