

Human PRCP protein (Recombinant) (His-Tag) (STJP019512)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid In Phosphate-Buffered Saline (pH 7.4) containing 30% Glycerol
Molecular Weight	54.3kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5547
Gene Symbol	PRCP
UniProt ID	PCP_HUMAN
Immunogen Region	22-496aa
Immunogen Sequence	LRPALRALGS LHLPTNPTSL PAVAKNYSVL YFQQKVDHFG FNTVKTFNQR YLVADKYWKK NGGSILFYTG NEGDIIWFCN NTGFMWDVAE ELKAMLVFAE HRYYGESLPF GDNSFKDSRH LNFLTSEQAL ADFAELIKHL KRTIPGAENQ PVIAIGGSYG GMLAAWFRMK YPHMVVGALA ASAPIWQFED LVPCGVFMKI VTDFRKS GP HCSESIHRSW DAINRLSN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human PRCP protein SDS-PAGE gel under reducing conditions stained with Coomassie blue