

Human DSPG3 protein (Recombinant) (His-Tag) (STJP019537)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Bioactivity
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	35.5kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1833
Gene Symbol	EPYC
UniProt ID	EPYC_HUMAN
Immunogen Region	20-322aa
Immunogen Sequence	APTL E SIN Y D D SET Y DAT L E D LD N LY N Y E NI I P V D K VE I E I A I TV M PS G N R E L L T PP P Q P E K A Q E E E E E E E E ST P RL I D G S S P Q E P E F T G V L G P HT N E D F P T C L L CT C I S T T V Y C D D H E L D A I P P L P K N T A Y F Y S R F N R I K K I N KN D F A S L . S D L KR I D L T S N L I I SE I D E D A F R K L P Q L R E L V L R D N K I R Q L P E L P T T L T F I D I S N N R L G R K G I K Q E A F K D M Y D L H H L Y L T D N

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human DSPG3 protein 3µg under reducing conditions with Coomassie blue stain