

Mouse CD99-L2 protein (Recombinant) (hIgG-His-Tag) (STJP019540)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	41.8kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	171486
Gene Symbol	Cd99l2
UniProt ID	C99L2_MOUSE
Immunogen Region	26-161aa
Immunogen Sequence	DTDGFNL EDALKETSSV KQRWDHFSTT TRRPVTTTRAP ANPAERWDHV ATTTTRRPGT TRAPSNPMEL DGFDELDALD DRNDLDGPKK PSAGEAGGWS DKDLEDIVEG GGYKPDKNKG GGGYGSNDDP GSGISTETG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Recombinant mouse CD99-L2 protein SDS-PAGE gel under reducing conditions stained with Coomassie blue