

Human IBSP/Sialoprotein II protein (Recombinant) (His-Tag) (STJP019542)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Bioactivity
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	34.3kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3381
Gene Symbol	IBSP
UniProt ID	SIAL_HUMAN
Immunogen Region	17-317aa
Immunogen Sequence	FSMKNLHRRV KIEDSEENG V FKYP RPYLY KHAYFYPHLK RFPVQGSSDS SEENGDDSS EEEEEETS N EGENNEESNE DEDSEAENTT LSATTLGYGE DATPGTGYTG LAAILPKKA GDITNKATKE KESDEEEEE EEGNENESE AEVDENEQGI NGTSTNSTEA ENGNSSGGD NGEEGEEESV TGANAEDTTE TGRQKGTSK TTTSPNGGFE PTTTPQVY

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human IBSP
Sialoprotein II protein under reducing conditions