

Human BAFF/TNFSF13B protein (Recombinant) (His-Tag) (STJP019551)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Bioactivity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) , 10% Glycerol
Molecular Weight	21.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10673
Gene Symbol	TNFSF13B
UniProt ID	TN13B_HUMAN
Immunogen Region	134-285aa
Immunogen Sequence	MAV QGPEETVTQD CLQLIADSET PTIQKGSYTF VPWLLSFKRG SALEEEKENKI LVKETGYFFI YGQVLYTDKT YAMGHLIQRK KVHVFGDELS LVTLFRCIQN MPETLPNNSC YSAGIAKLEE GDELQLAIPR ENAIQSLDGD VTFFGALKLL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human BAFF TNFSF13B protein 3ug under reducing conditions Coomassie blue stain



Human BAFF protein ELISA binding assay with BAFFR receptor