

Human TMIGD2/CD28H protein (Recombinant) (hIgG-Tag) (STJP019554)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Bioactivity
Host / Source	HEK293 Cells

PRODUCT PROPERTIES



Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	40.1kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	126259
Gene Symbol	TMIGD2
UniProt ID	TMIG2_HUMAN
Immunogen Region	23-150aa
Immunogen Sequence	LSVQQGPNLL QVRQGSQATL VCQVDQATAW ERLRVKWKTD GAILCQPYIT NGSLSLGVCG PQGRLSWQAP SHLTQLQDPV SLNHSGAYVC WAAVEIPELE EAEGNITRLF VDPDDPTQNR NRIASFPG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel electrophoresis of recombinant human TMIGD2 CD28H protein under reducing conditions

Human TMIGD2 CD28H protein binding to HHLA2 biological activity assay