

Human GFR alpha-1 protein (Recombinant) (hlgG-Tag) (STJP019564)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Bioactivity
Host / Source	HEK293 Cells

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	70.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2674
Gene Symbol	GFRA1
UniProt ID	GFRA1_HUMAN
Immunogen Region	25-424aa
Immunogen Sequence	DRLDCKVQASD QCLKEQSCST KYRTLRCQVA GKETNFSLAS GLEAKDECERS AMEALKQKSL YNCRCCKRGMK KEKNCLRIYW SMYQSLQGND LLEDSPYEPV NSRLSDIFRV VPFISVEHIP KGNNCLDAK ACNLDDICKK YRSAYITPCT TSVSNDVCNR RKCHKALRQF FDKVPAKHSY GMLFCSCRD IACTERRRQTI VPVCSYEERE KPNCLNLQDS CKTNYICR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human GFR alpha-1 protein under reducing conditions stained with Coomassie blue