

## Human Fc gamma RIIA/CD32a (R167) protein (Recombinant) (His-Tag) (STJP019569)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Bioactivity
Host / Source	HEK293 Cells

### PRODUCT PROPERTIES



Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	21.7kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2212</a>
Gene Symbol	<a href="#">FCGR2A</a>
UniProt ID	<a href="#">FCG2A_HUMAN</a>
Immunogen Region	30-218aa
Immunogen Sequence	SADSQAAAPP KAVLKLEPPW INVLQEDSVT LTCQGARSPE SDSIQWFHNG NLIPTHTQPS YRFKANNNDS GEYTCQTGQT SLSDPVHLTV LSEWLVLTQP HLEFQEGETI MLRCHSWKDK PLVKVTFEQN GKSKQKFSRLD PTFISIPQANH SHSGDYHCTG NIGYTLFSSK PVTITVQVPS MGSSSPMGI

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human Fc gamma RIIA CD32a R167 protein under reducing conditions  
Fc gamma RIIA CD32a R167 protein functional ELISA binding assay with human IgG