

Human BID protein (Recombinant) (No-Tag) (STJP019572)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol
Molecular Weight	21.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	637
Gene Symbol	BID
UniProt ID	BID_HUMAN
Immunogen Region	1-195aa
Immunogen Sequence	MDCEVNNQSS LRDECITNLL VFGFLQSCSD NSFRRDLAL GHLPVLAPQ WEGYDELQTD GNRSSHSRLG RIEADSESQE DIIRNIARHL AQVGDSMDRS IPPGLVNGLA LQLRNTSRSE EDNRNRLATA LEQLLQAYPR DMEKEKTMVL LALLAKKVA SHTPSLLRDV FHTTVNFINQ NLRITYVRS LA RINGMD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human BID protein 3µg under reducing conditions with Coomassie blue stain