

Human BAFF/TNFSF13B protein (Recombinant) (His-Tag) (STJP019576)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) , 10% Glycerol
Molecular Weight	21.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10673
Gene Symbol	TNFSF13B
UniProt ID	TN13B_HUMAN
Immunogen Region	134-285aa
Immunogen Sequence	MAV QGPEETVTQD CLQLIADSET PTIQKGSYTF VPWLLSFKRG SALEEEKENKI LVKETGYFFI YGQVLYTDKT YAMGHLIQRK KVHVFGDELS LVTLFRCIQN MPETLPNNNSC YSAGIAKLEE GDELQLAIPR ENAIQISLDGD VTFFGALKLL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing 3µg recombinant human BAFF TNFSF13B protein under reducing conditions Coomassie blue stain