

Human BMP-7 (monomer) protein (Recombinant) (His-Tag) (STJP019579)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 10mM Sodium Citrate buffer (pH 3.5) containing 10% Glycerol
Molecular Weight	16.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	655
Gene Symbol	BMP7
UniProt ID	BMP7_HUMAN
Immunogen Region	293-431aa
Immunogen Sequence	MSTGSKQRSQ NRSKTPKNQE ALRMANVAEN SSSDQRQACK KHELYVSFRD LGWQDWIIAP EGYAAYYCEG ECAFLPLNSYM NATNHAIVQT LVHFINPETV PKPCCAPTQL NAISVLYFDD SSNVILKKYR NMVVRACGCH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human BMP-7 monomer protein under reducing conditions stained with Coomassie blue