

Human BNP protein (Recombinant) (His-Tag) (STJP019580)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) 40% Glycerol, 2mM DTT, 0.1mM PMSF, 1mM EDTA
Molecular Weight	14 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4879
Gene Symbol	NPPB
UniProt ID	ANFB_HUMAN
Immunogen Region	27-134aa
Immunogen Sequence	HPLGSPGSA SDLETSGLQE QRNHLQGKLS ELQVEQTSLE PLQESPRPTG VWKSREVATE GIRGHRKMVL YTLRAPRSPK MVQGSGCFGR KMDRISSSSG LGCKVLRRH

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human BNP protein 3µg under reducing conditions Coomassie blue stain