

Human CEBP alpha/CEBPA protein (Recombinant) (His-Tag) (STJP019584)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 0.1M NaCl, 5mM beta-Mercaptoethanol.
Molecular Weight	14.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1050
Gene Symbol	CEBPA
UniProt ID	CEBPA_HUMAN
Immunogen Region	270-358aa
Immunogen Sequence	GAG KAKKSVDKNS NEYRVRERIN NIAVRKSRDK AKQRNVETQQ KVLELTSDND RLRKRVEQLS RELDTLRGIF RQLPESSLVK AMGNCA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human CEBPA protein 3µg under reducing conditions Coomassie blue stain