

Human HABP1/C1QBP protein (Recombinant) (No-Tag) (STJP019594)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 1mM DTT, 20% Glycerol
Molecular Weight	23.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	708
Gene Symbol	C1QBP
UniProt ID	C1QBP_HUMAN
Immunogen Region	74-282aa
Immunogen Sequence	MLHTDGDKAF VDFLSDEIKE ERKIQKHKTL PKMSGGWELE LNGTEAKLVR KVAGEKITVT FNINNSIPPT FDGEEEPSQG QKVVEEQPEL TSTPNFVVEV IKNDDGKKAL VLDCHYPEDE VGQEDEAESD IFSIREVSFQ STGESEWKDT NYTLNTDSL D WALYDHLMDF LADRGVDNTF ADELVELSTA LEHQEYITFL EDLKSFVKSQ

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human HABP1 C1QBP protein under reducing conditions