

## Human CRADD protein (Recombinant) (His-Tag) (STJP019595)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol 1mM DTT
Molecular Weight	24.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">8738</a>
Gene Symbol	<a href="#">CRADD</a>
UniProt ID	<a href="#">CRADD_HUMAN</a>
Immunogen Region	1-199aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MEARDKQVLR SLRLELGAEV LVEGLVLQYL YQEGILTENH IQEINAQTTG LRKTMMLLDI LPSRGPKAFD TFLDSLQEFV WVREKLLKAR EEAMTDLPAQ DRLTGIPSHI LNSSPSDRQI NQLAQRGLGPE WPEMVLSLGL SQTDIYRCKA NHPHNVSQV VEA FIRWRQR FGKQATFQSL HNGLRAVEVD PSLLLHMLE

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human CRADD protein under reducing conditions