

Human GRO gamma/CXCL3 protein (Recombinant) (His-Tag) (STJP019602)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 20% Glycerol
Molecular Weight	10.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2921
Gene Symbol	CXCL3
UniProt ID	CXCL3_HUMAN
Immunogen Region	35-107aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MASVVTCLRC QCLQTLQGIH LKNIQSVNVR SPGPHCAQTE VIATLKNGKK ACLNPAAPMV QKIIEKILNK GSTN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human GRO gamma CXCL3 protein 3 µg reducing conditions
Coomassie blue