

Human Cyclophilin B/PPIB protein (Recombinant) (No-Tag) (STJP019603)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20mM NaCl, 0.5mM DTT, 10% Glycerol
Molecular Weight	21.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5479
Gene Symbol	PPIB
UniProt ID	PPIB_HUMAN
Immunogen Region	26-216aa
Immunogen Sequence	MLLPGPSAAD EKKKGPKVTV KVFYDLRIGD EDVGRVIFGL FGKTVPKTVD NFVALATGEK GFGYKNSKFH RVIKDFMIQG GDFTRGDGTG GKSİYGERFP DENFKLKHYG PGWVSMANAG KDTNGSQFFI TTVKTAWLDG KHVVFQKVL E GMEVVRKVES TKTDSRDKPL KDVIADCGK IVEKPFIAIA KE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Cyclophilin B PPIB protein under reducing conditions