


E.coli DnaK (508-638aa) protein (Recombinant) (No-Tag) (STJP019613)

GENERAL INFORMATION		
	Product Type	Proteins
	Applications	SDS-PAGE
	Host / Source	E.coli
PRODUCT PROPERTIES		
	Concentration	1 mg/mL
	Formulation	Liquid in 25mM Tris-HCl buffer (pH 7.5) containing 100mM NaCl, 10% Glycerol
	Molecular Weight	14.6 kDa
	Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Immunogen Region	508-638aa
	Immunogen Sequence	MNEDEIQKMV RDAEANA EAD RKFEELVQTR NQGDHLLHST RKQVEEAGDK LPADDKTAIE SALTALETAL KGEDKAAIEA KMQELAQVSQ KLMEIAQQQH AQQQTAGADA SANNAKDDDV VDAEFEEVKD KK
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant E.coli DnaK 508-638aa protein under reducing conditions