

Human Calcipressin-1/DSCR1 protein (Recombinant) (No-Tag) (STJP019614)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 0.1M NaCl, 1mM DTT, 10% Glycerol
Molecular Weight	13.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1827
Gene Symbol	RCAN1
UniProt ID	RCAN1_HUMAN
Immunogen Region	1-117aa
Immunogen Sequence	MKLYFAQTLH IGSSHLAPPN PDKQFLISPP ASPPVGWKQV EDATPVINYD LLYAISKLGP GEKYELHAAT DTTSPVVVHV CESDQEKEEE EEMERMRRPK PKIIQTRRPE YTIPIHLS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Calcipressin-1 DSCR1 protein 3µg under reducing conditions Coomassie blue stain