

Human EBAG9 protein (Recombinant) (His-Tag) (STJP019615)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 2mM EDTA, 10% Glycerol
Molecular Weight	23.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9166
Gene Symbol	EBAG9
UniProt ID	RCAS1_HUMAN
Immunogen Region	28-213aa
Immunogen Sequence	RSGRGRKLS GDQITLPTTV DYSSVPKQTD VEEWTSWDED APTSVKIEGG NGNVATQQNS LEQLEPDYFK DMTPTIKKTQ KIVIKKREPL NFGIPDGSTG FSSRLAATQD LPFIHQSSSEL GDLDTWQENT NAWEEEEEDAA WQAEVLRQQ KLADREKRAA EQQRKKMEKE AQRLMKKEQN KIGVKLS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing 3µg recombinant human EBAG9 protein under reducing conditions stained with Coomassie blue