

## Human FABP7/B-FABP protein (Recombinant) (No-Tag) (STJP019619)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 25mM Tris-HCl buffer (pH 7.5) containing 1mM DTT, 2mM EDTA, 10% Glycerol
Molecular Weight	14 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2173</a>
Gene Symbol	<a href="#">FABP7</a>
UniProt ID	<a href="#">FABP7_HUMAN</a>
Immunogen Region	1-132aa
Immunogen Sequence	MVEAFCATWK LTNSQNFDEY MKALGVGFAT RQVGNVTKPT VIISQEGDKV VIRTSTFKN TEISFQLGEE FDETTADDRN CKSVVSLDGD KLVHIQKWDG KETNFVREIK DGKMMVMTLTF GDVVAVRHYE KA

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human FABP7 B-FABP protein STJP019619 under reducing conditions