

## Human FABP4/A-FABP protein (Recombinant) (No-Tag) (STJP019620)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 2mM EDTA, 10% Glycerol
Molecular Weight	14.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2167</a>
Gene Symbol	<a href="#">FABP4</a>
UniProt ID	<a href="#">FABP4_HUMAN</a>
Immunogen Region	1-132aa
Immunogen Sequence	MCDAFVGTWK LVSSNFDDY MKEVGVGFAT RKVAGMAKPN MIISVNGDVI TIKSESTFKN TEISFILGQE FDEVTADDRK VKSTITLDGG VLVHVQKWDG KSTTIKRE DDKLVVECVM KGVSTRVYE RA

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human FABP4 A-FABP protein 3µg under reducing conditions with Coomassie blue stain