

Human FABP5/E-FABP protein (Recombinant) (No-Tag) (STJP019623)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 20% Glycerol
Molecular Weight	15.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2171
Gene Symbol	FABP5
UniProt ID	FABP5_HUMAN
Immunogen Region	1-135aa
Immunogen Sequence	MATVQQLEGR WRLVDSKGF D EYMKELGVGI ALRKMGMAMAK PDCIITCDGK NLTIKTESTL KTTQFSCTLG EKFEETTADG RKTQTVCNFT DGALVQHGEW DGKESTITRK LKDGKLVVEC VMNNVTCTRI YEKVE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing 3µg recombinant human FABP5 E-FABP protein under reducing conditions
 Coomassie blue stain