

Human FTL protein (Recombinant) (No-Tag) (STJP019628)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5)
Molecular Weight	20kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2512
Gene Symbol	FTL
UniProt ID	FRIL_HUMAN
Immunogen Region	1-175aa
Immunogen Sequence	MSSQIRQNYSDVEAAVNSLVNLYLQASYTYLSLGFYFDRDDVALEGVSHFFRELAAEEKREGYERLLKMQNQRGGRALFQDIKKPAEDEWGKTPDAMKAA MALEKKLNQA LLDLHALGSA RTDPHLCDFL ETHFLDEEVK LIKKMGDHLT NLHRLGGPEA GLGEYLFERL TLKHD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human FTL protein 3µg under reducing conditions stained with Coomassie blue