

Human GMF-gamma protein (Recombinant) (No-Tag) (STJP019637)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 1mM EDTA, 10% Glycerol
Molecular Weight	16.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9535
Gene Symbol	GMFG
UniProt ID	GMFG_HUMAN
Immunogen Region	1-142aa
Immunogen Sequence	MSDSLVVCEV DPTELTKLRK FRFRKETDNA AIIMKVDKDR QMVLLEEEFQ NISPEELKME LPERQPRFVV YSYKYVHDDG RVSYP LCFIF SSPVGCKPEQ QMMYAGSKNR LVQTAELTKV FEIRTTDDL T EAWLQEKLSF FR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human GMF-gamma protein 3 μ g reducing conditions
Coomassie blue stain