

e.coli GroES protein (Recombinant) (No-Tag) (STJP019645)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 25mM Tris-HCl buffer (pH 7.5) containing 100mM NaCl, 5mM DTT, 10%Glycerol
Molecular Weight	10.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	1-97aa
Immunogen Sequence	MNIRPLHDRV IVKRKEVETK SAGGIVLTS AAKSTRGEV LAVNGRILE NGEVKPLDVK VGDIVFNDG YGVKSEKIDN EEVLIMSED ILAIVEA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant E. coli GroES protein under reducing conditions