

Human HIP/ST13 protein (Recombinant) (No-Tag) (STJP019651)

GENERAL INFORMATION

| | |
|---------------|----------|
| Product Type | Proteins |
| Applications | SDS-PAGE |
| Host / Source | E.coli |

PRODUCT PROPERTIES


| | |
|---------------------|---|
| Concentration | 1 mg/mL |
| Formulation | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 0.1M NaCl, 10% Glycerol |
| Molecular Weight | 41.3 kDa |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|--------------------|--|
| Gene ID | 6767 |
| Gene Symbol | ST13 |
| UniProt ID | F10A1_HUMAN |
| Immunogen Region | 1-369aa |
| Immunogen Sequence | MDPRKVNELR AFVKMCKQDP SVLHTEEMRF LEWVESMGG KVPPATQKAK SEENTKEEKP DSKKVEEDLK ADEPSSEESD LEIDKEGVIE PDTDAPQEMG DENAEITEEM MDQANDKKVA AIEALNDGEL QKAILDFTDA IKLNPRILAIL YAKRASVFKV LQKPNAAIRD CDRAIEINPD SAQPYKWRGK AHRLLGHWEE AAHDLALACK LDYDEASAM LKEVQPRA |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human HIP ST13 protein 3µg Coomassie blue stain