

Human HSBP1 protein (Recombinant) (No-Tag) (STJP019654)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing, 50mM NaCl, 20% Glycerol
Molecular Weight	8.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3281
Gene Symbol	HSBP1
UniProt ID	HSBP1_HUMAN
Immunogen Region	1-76aa
Immunogen Sequence	MAETDPKTVQ DLTSVVQTLL QQMQDKFQTM SDQIIIGRIDD MSSRIDDLK NIADLMTQAG VEELESENKI PATQKS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human HSBP1 protein under reducing conditions stained with Coomassie blue