

Human HSP90 alpha protein (Recombinant) (His-Tag) (STJP019656)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid In Phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	86.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3320
Gene Symbol	HSP90AA1
UniProt ID	HS90A_HUMAN
Immunogen Region	1-732aa
Immunogen Sequence	MPEETQTQDQ PMEEEEVETF AFQAEIAQLM SLIINTFYNS KEIFLRELIS NSSDALDKIR YESLTDPSKL DSGKELHINL IPNKQDRTLIT IVDTGIGMTK ADLINNLGTI AKSGTKAFME ALQAGADISM IGQFGVGFYS AYLVAEKVTV ITKHNDDEQY AWESSAGGSF TVRTDTGEPM GRGTKVILHL KEDQTEYLEE RRIKEI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human HSP90 alpha protein 3 µg under reducing conditions with Coomassie blue staining