

Human HSP40/DNAJB1 protein (Recombinant) (His-Tag) (STJP019662)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	40.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3337
Gene Symbol	DNAJB1
UniProt ID	DNJB1_HUMAN
Immunogen Region	1-340aa
Immunogen Sequence	MGKDYYQTLG LARGASDEEI KRAYRRQALR YHPDKNKEPG AEEKFKEIAE AYDVLSDPRK REIFDRYGEELGLKGSGPSGG SGGGANGTSF SYTFHGDPHA MFAEFFGGRN PFDTFFGQRN GEEGMDIDDPFSGFPMGMGG FTNVNFGRRS SAQEPARKKQ DPPVTHDLRV SLEEIYSGCT KMKISHKRL NPDGKSIRNE DKILTI

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human HSP40/DNAJB1 protein under reducing conditions stained with Coomassie blue